

# Resource Summary Report

Generated by [kravitz2](#) on Jul 30, 2026

## Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size

RRID:SCR\_020040

Type: Tool

### Proper Citation

Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size (RRID:SCR\_020040)

### Resource Information

**URL:** [https://www.horiba.com/en\\_en/products/detail/action/show/Product/camsizer-x2-1133/](https://www.horiba.com/en_en/products/detail/action/show/Product/camsizer-x2-1133/)

**Proper Citation:** Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size (RRID:SCR\_020040)

**Description:** Dynamic image analysis system that provides particle size and shape distribution analysis needed for quality control of fine powders and suspensions, measuring 0.8  $\mu\text{m}$  to 8 mm.

**Resource Type:** instrument resource

**Keywords:** Horiba, Dynamic Image Analysis System Particle Size, Instrument Equipment, USEDit,

**Availability:** Commercially available

**Resource Name:** Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size

**Resource ID:** SCR\_020040

**Alternate IDs:** Model\_Number\_Camsizer X2

**Record Creation Time:** 20220129T080348+0000

**Record Last Update:** 20260725T190920+0000

### Ratings and Alerts

No rating or validation information has been found for Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size.

No alerts have been found for Horiba: Camsizer X2 Dynamic Image Analysis System Particle Size.

---

## Data and Source Information

**Source:** [SciCrunch Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.