

Resource Summary Report

Generated by [nidm-terms](#) on Jul 28, 2026

Thermo Fisher: Applied Biosystems: ABI Voyager DE Pro MALDI-Mass Spectrometer

RRID:SCR_019340

Type: Tool

Proper Citation

Thermo Fisher: Applied Biosystems: ABI Voyager DE Pro MALDI-Mass Spectrometer (RRID:SCR_019340)

Resource Information

URL: <https://www.labwrench.com/equipment/8044/applied-biosystems-voyager-de-trade-pro>

Proper Citation: Thermo Fisher: Applied Biosystems: ABI Voyager DE Pro MALDI-Mass Spectrometer (RRID:SCR_019340)

Description: This mass spectrometer can be used for structural analysis of organic and inorganic compounds, protein identification studies, synthetic peptide validation experiments, small molecule and polymer analysis, and post source decay experiments.

Resource Type: instrument resource

Keywords: ABI, Applied Biosystems, Mass Spectrometer, Matrix Assisted Laser Desorption Ionization Time of Flight Instrument, Instrument Equipment, USEDit

Availability: Commercially available

Resource Name: Thermo Fisher: Applied Biosystems: ABI Voyager DE Pro MALDI-Mass Spectrometer

Resource ID: SCR_019340

Alternate IDs: Model_Number_Voyager

Record Creation Time: 20220129T080344+0000

Record Last Update: 20260725T190908+0000

Ratings and Alerts

No rating or validation information has been found for Thermo Fisher: Applied Biosystems: ABI Voyager DE Pro MALDI-Mass Spectrometer.

No alerts have been found for Thermo Fisher: Applied Biosystems: ABI Voyager DE Pro MALDI-Mass Spectrometer.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We found 1 mentions in open access literature.

Listed below are recent publications. The full list is available at [nidm-terms](#) .

Cann P, et al. (2022) Variation of ewe olfactory secretome during a ram effect. Frontiers in veterinary science, 9, 1033412.