

# Resource Summary Report

Generated by [PRECISE-TBI](#) on Jul 31, 2026

## Sage Science: Pippin Prep

RRID:SCR\_020403

Type: Tool

---

### Proper Citation

Sage Science: Pippin Prep (RRID:SCR\_020403)

---

### Resource Information

**URL:** <https://sagescience.com/products/pippin-prep/>

**Proper Citation:** Sage Science: Pippin Prep (RRID:SCR\_020403)

**Description:** Facilitates library construction system for NGS platforms, and is recommended by Illumina and Ion Torrent for certain workflows. The platform features the ability to collect narrow and even fragment distributions, as well as the flexibility to collect wide ranges of fragments with minimal effort. Target sizes or ranges of sizes are entered in software, and fractions are collected in buffer. Up to 5 samples per gel cassette may be run, with no possibility of cross contamination.

**Resource Type:** instrument resource

**Keywords:** Library construction, NGS platform, instrument, equipment, USEDit

**Availability:** Commercially available

**Resource Name:** Sage Science: Pippin Prep

**Resource ID:** SCR\_020403

**Alternate IDs:** Model\_Number\_Pippin Prep

**Record Creation Time:** 20220129T080350+0000

**Record Last Update:** 20260725T190928+0000

---

### Ratings and Alerts

No rating or validation information has been found for Sage Science: Pippin Prep.

No alerts have been found for Sage Science: Pippin Prep.

---

### Data and Source Information

**Source:** [SciCrunch Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.