

Resource Summary Report

Generated by T1D on Jul 31, 2026

Beckman Coulter: ACT10 Cell differential counter

RRID:SCR_019621

Type: Tool

Proper Citation

Beckman Coulter: ACT10 Cell differential counter (RRID:SCR_019621)

Resource Information

URL: <https://analyticalinstrument.com/index.php/lab-equipment/beckman-coulter-act-10-6706319-detail>

Proper Citation: Beckman Coulter: ACT10 Cell differential counter (RRID:SCR_019621)

Description: AcT 10 is an and hematology analyzer with whole blood and pre-dilute capabilities from Beckman/Coulter. Beckman/Coulter AcT 10 automated analyzer is intended for testing in any smaller lab facility and low-volume office labs. Ten parameters including WBC, RBC, Hgb, Hct, MCV, MCH, MCHC, Plt, Lymph # and Lymph % can be measured using this Beckman/Coulter AcT 10 analyzer. For abnormal samples, partial differential option is available. Sample volume required for the whole blood mode is 12 ?L and 20 ?L in pre-dilute mode.

Resource Type: instrument resource

Keywords: Beckman Coulter, Cell Differential Counter, Instrument Equipment, USEDit

Availability: Commercially available

Resource Name: Beckman Coulter: ACT10 Cell differential counter

Resource ID: SCR_019621

Alternate IDs: Model_Number_ACT10

Record Creation Time: 20220129T080346+0000

Record Last Update: 20260725T190913+0000

Ratings and Alerts

No rating or validation information has been found for Beckman Coulter: ACT10 Cell differential counter.

No alerts have been found for Beckman Coulter: ACT10 Cell differential counter.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.