

# Resource Summary Report

Generated by T1D on Aug 26, 2026

## Bruker: timsTOF Ultra 2 mass spectrometer

RRID:SCR\_026541

Type: Tool

### Proper Citation

Bruker: timsTOF Ultra 2 mass spectrometer (RRID:SCR\_026541)

### Resource Information

**URL:** <https://www.bruker.com/content/bruker/int/en/products-and-solutions/mass-spectrometry/timstof/timstof-ultra.html>

**Proper Citation:** Bruker: timsTOF Ultra 2 mass spectrometer (RRID:SCR\_026541)

**Description:** Spectrometer timsTOF Ultra 2 is paired with CaptiveSpray Ultra 2 ion source for high sensitivity to work with challenging samples and detect low-abundance peptides. Used to analyze even the smallest amounts and uncover previously unseen peptides.

**Synonyms:** , Bruker timsTOF Ultra 2, timsTOF Ultra 2 mass spectrometer

**Resource Type:** instrument resource

**Keywords:** CaptiveSpray Ultra 2 ion source, high sensitivity, challenging samples, detect low-abundance peptides, unseen peptides,

**Specification URL:** [https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/main/PDF/SCR\\_026541.pdf](https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/main/PDF/SCR_026541.pdf)

**Resource Name:** Bruker: timsTOF Ultra 2 mass spectrometer

**Resource ID:** SCR\_026541

**Alternate IDs:** Model\_Number\_timsTOF\_Ultra\_2

**Record Creation Time:** 20250307T060316+0000

**Record Last Update:** 20260815T183242+0000

### Ratings and Alerts

No rating or validation information has been found for Bruker: timsTOF Ultra 2 mass spectrometer.

No alerts have been found for Bruker: timsTOF Ultra 2 mass spectrometer.

---

## Data and Source Information

**Source:** [SciCrunch Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.