

Resource Summary Report

Generated by T1D on Aug 12, 2026

Planmeca: Emerald Intraoral scanner

RRID:SCR_018404

Type: Tool

Proper Citation

Planmeca: Emerald Intraoral scanner (RRID:SCR_018404)

Resource Information

URL: <https://www.planmeca.com/cadcam/dental-scanning/planmeca-emerald/>

Proper Citation: Planmeca: Emerald Intraoral scanner (RRID:SCR_018404)

Description: Intraoral scanner for taking digital impressions.

Synonyms: Planmeca Emerald, Planmeca Emerald TM Intraoral scanner

Resource Type: instrument resource

Keywords: Intraoral scanner, digital impression, scanner, instrument, equipment, USEDit

Availability: Restricted

Resource Name: Planmeca: Emerald Intraoral scanner

Resource ID: SCR_018404

Alternate IDs: Model_Number_Emerald

Alternate URLs: <https://www.planmeca.com/globalassets/planmeca/1-planmeca-usa-modernization/cadcam/dental-scanning/emerald/manuals/emerald-manual-15698200e.pdf>

Record Creation Time: 20220129T080340+0000

Record Last Update: 20260808T190105+0000

Ratings and Alerts

No rating or validation information has been found for Planmeca: Emerald Intraoral scanner.

No alerts have been found for Planmeca: Emerald Intraoral scanner.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.