

Resource Summary Report

Generated by [kravitz2](#) on Sep 5, 2026

Donkey anti-Goat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 550

RRID:AB_2556667

Type: Antibody

Proper Citation

Thermo Fisher Scientific Cat# SA5-10087, RRID:AB_2556667

Antibody Information

URL: http://antibodyregistry.org/AB_2556667

Proper Citation: Thermo Fisher Scientific Cat# SA5-10087, RRID:AB_2556667

Target Antigen: Goat IgG (H+L)

Host Organism: donkey

Clonality: polyclonal secondary

Comments: Applications: Flow, ICC/IF, IHC, IP, WB

Antibody Name: Donkey anti-Goat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 550

Description: This polyclonal secondary targets Goat IgG (H+L)

Target Organism: goat

Antibody ID: AB_2556667

Vendor: Thermo Fisher Scientific

Catalog Number: SA5-10087

Record Creation Time: 20231110T035234+0000

Record Last Update: 20260829T152044+0000

Ratings and Alerts

No rating or validation information has been found for Donkey anti-Goat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 550.

No alerts have been found for Donkey anti-Goat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 550.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 6 mentions in open access literature.

Listed below are recent publications. The full list is available at [kravitz2](#) .

Andreska T, et al. (2023) DRD1 signaling modulates TrkB turnover and BDNF sensitivity in direct pathway striatal medium spiny neurons. Cell reports, 42(6), 112575.

McLeod VM, et al. (2022) Mapping Motor Neuron Vulnerability in the Neuraxis of Male SOD1G93A Mice Reveals Widespread Loss of Androgen Receptor Occurring Early in Spinal Motor Neurons. Frontiers in endocrinology, 13, 808479.

Tomas D, et al. (2021) Dissociation of disease onset, progression and sex differences from androgen receptor levels in a mouse model of amyotrophic lateral sclerosis. Scientific reports, 11(1), 9255.

McLeod VM, et al. (2019) Androgen receptor antagonism accelerates disease onset in the SOD1G93A mouse model of amyotrophic lateral sclerosis. British journal of pharmacology, 176(13), 2111.

Benz F, et al. (2019) Low wnt/β-catenin signaling determines leaky vessels in the subfornical organ and affects water homeostasis in mice. eLife, 8.

Li R, et al. (2019) Generation of Blastocyst-like Structures from Mouse Embryonic and Adult Cell Cultures. Cell, 179(3), 687.