

Resource Summary Report

Generated by [kravitz2](#) on Jul 30, 2026

Thermo Fisher: Dionex PC10 Postcolumn Pneumatic Delivery

RRID:SCR_019841

Type: Tool

Proper Citation

Thermo Fisher: Dionex PC10 Postcolumn Pneumatic Delivery (RRID:SCR_019841)

Resource Information

URL: <https://www.thermofisher.com/order/catalog/product/053591>

Proper Citation: Thermo Fisher: Dionex PC10 Postcolumn Pneumatic Delivery (RRID:SCR_019841)

Description: Reaction system that ensures consistent, pulse free reagent delivery for applications with high background absorbances, such as trace level transition metal analysis using post column reaction with 4 2 pyridylazoresorcinol PAR. Provides reagent with less pressure pulsation by applying constant pressure, eliminating need for pump. Assures optimal flow by applying constant pressure, eliminating need for pump.

Synonyms: PC10 Postcolumn Pneumatic Delivery System

Resource Type: instrument resource

Keywords: Thermo Fisher, Thermo Fisher, Ion Chromatography Component, Instrument, Equipment, USEDit,

Availability: Commercially available

Specification URL: https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/main/PDF/SCR_019841.pdf

Resource Name: Thermo Fisher: Dionex PC10 Postcolumn Pneumatic Delivery

Resource ID: SCR_019841

Alternate IDs: Model_Number_Dionex_PC10

Record Creation Time: 20220129T080347+0000

Record Last Update: 20260725T190919+0000

Ratings and Alerts

No rating or validation information has been found for Thermo Fisher: Dionex PC10 Postcolumn Pneumatic Delivery.

No alerts have been found for Thermo Fisher: Dionex PC10 Postcolumn Pneumatic Delivery.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.