

Resource Summary Report

Generated by [kravitz2](#) on Jul 29, 2026

Horiba: Fluorolog-3 Spectrofluorometer

RRID:SCR_020061

Type: Tool

Proper Citation

Horiba: Fluorolog-3 Spectrofluorometer (RRID:SCR_020061)

Resource Information

URL: https://www.horiba.com/en_en/products/detail/action/show/Product/fluorolog-1586/

Proper Citation: Horiba: Fluorolog-3 Spectrofluorometer (RRID:SCR_020061)

Description: Fluorolog-3 is a , modular system which allows the researcher to interchange a wide range of accessories to correspond precisely with the characteristics of a given sample. From analysis of steady-state or molecular dynamics to IR probes, the Fluorolog-3 comes equipped with a wide range and limitless configuration of accessories to enhance the accuracy and speed of it's performance. The Fluorolog-3 series enables you to customize your spectrofluorometer's performance to meet your research needs. It can also perform (TCSPC) Time Correlated Single Photon Counting lifetime measurements like the FluoroMax-4 and QM800.

Resource Type: instrument resource

Keywords: Horiba, Spectrofluorometer, Instrument Equipment, USEdit

Availability: Commercially available

Resource Name: Horiba: Fluorolog-3 Spectrofluorometer

Resource ID: SCR_020061

Alternate IDs: Model_Number_Fluorolog-3

Record Creation Time: 20220129T080348+0000

Record Last Update: 20260725T190923+0000

Ratings and Alerts

No rating or validation information has been found for Horiba: Fluorolog-3 Spectrofluorometer.

No alerts have been found for Horiba: Fluorolog-3 Spectrofluorometer.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We found 1 mentions in open access literature.

Listed below are recent publications. The full list is available at [kravitz2](#) .

Seraphim TV, et al. (2022) Assembly principles of the human R2TP chaperone complex reveal the presence of R2T and R2P complexes. Structure (London, England : 1993), 30(1), 156.