

Resource Summary Report

Generated by [kravitz2](#) on Aug 5, 2026

Bruker: timsTOF fleX with MALDI-2

RRID:SCR_023615

Type: Tool

Proper Citation

Bruker: timsTOF fleX with MALDI-2 (RRID:SCR_023615)

Resource Information

URL: <https://www.bruker.com/en/products-and-solutions/mass-spectrometry/timstof/timstof-flex-maldi-2.html>

Proper Citation: Bruker: timsTOF fleX with MALDI-2 (RRID:SCR_023615)

Description: Instrument for probing molecular content of sample with MALDI-2. The timsTOF fleX is fully functional high speed, high sensitivity ESI instrument for all X-Omics analyses, combined with high spatial resolution Matrix-Assisted Laser Desorption Ionization (MALDI) source and stage specifically designed for resolving molecular distributions and bringing spatial dimension to OMICS analyses. MALDI-2 uses laser post-ionization to enhance and enrich the MALDI experiment providing access to chemical classes typically opaque to MALDI with enhanced sensitivities.

Resource Type: instrument resource

Keywords: Thermal Ionization Mass Spectrometry, probing molecular content, X-Omics analyses, instrument, equipment, USEDit

Availability: Commercially available

Specification URL: https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/refs/heads/main/PDF/SCR_023615.pdf

Resource Name: Bruker: timsTOF fleX with MALDI-2

Resource ID: SCR_023615

Record Creation Time: 20230527T050216+0000

Record Last Update: 20260801T190821+0000

Ratings and Alerts

No rating or validation information has been found for Bruker: timsTOF fleX with MALDI-2.

No alerts have been found for Bruker: timsTOF fleX with MALDI-2.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.