

Resource Summary Report

Generated by [nidm-terms](#) on Aug 10, 2026

Goat anti-Rabbit IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor™ Plus 488

RRID:AB_2633280

Type: Antibody

Proper Citation

(Thermo Fisher Scientific Cat# A32731, RRID:AB_2633280)

Antibody Information

URL: http://antibodyregistry.org/AB_2633280

Proper Citation: (Thermo Fisher Scientific Cat# A32731, RRID:AB_2633280)

Target Antigen: Rabbit IgG (H+L)

Host Organism: goat

Clonality: polyclonal secondary

Comments: Applications: ICC/IF, WB

Antibody Name: Goat anti-Rabbit IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor™ Plus 488

Description: This polyclonal secondary targets Rabbit IgG (H+L)

Target Organism: rabbit

Antibody ID: AB_2633280

Vendor: Thermo Fisher Scientific

Catalog Number: A32731

Record Creation Time: 20251216T074310+0000

Record Last Update: 20260225T231712+0000

Ratings and Alerts

No rating or validation information has been found for Goat anti-Rabbit IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor™ Plus 488.

No alerts have been found for Goat anti-Rabbit IgG (H+L) Highly Cross-Adsorbed Secondary Antibody, Alexa Fluor™ Plus 488.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 386 mentions in open access literature.

Listed below are recent publications. The full list is available at [nidm-terms](#) .

Tsutsumi H, et al. (2026) Single-nucleus transcriptional and chromatin accessibility analyses of maturing mouse Achilles tendon uncover the molecular landscape of tendon stem/progenitor cells. *eLife*, 14.

Fischer S, et al. (2026) Loss of ADAMTS9 disrupts ciliogenesis and collagen homeostasis resulting in Nephronophthisis-like polycystic kidneys. *bioRxiv : the preprint server for biology*.

de Vries DK, et al. (2026) Generation of the human induced pluripotent stem cell (hiPSC) line AUMCi015-A from a heterozygous carrier of the LMNA p.Q493X rare pathogenic variant. *Stem cell research*, 92, 103962.

Alunni A, et al. (2026) A549 and Ci-hAELVi cell lines coculture as a new human alveolar epithelium model for air-liquid interface exposure. *Toxicology in vitro : an international journal published in association with BIBRA*, 112, 106200.

Karnay A, et al. (2026) Transcription and cohesin direct domain boundary spatial positioning and are linked to Friedreich's ataxia. *Molecular cell*, 86(10), 1911.

Fatima M, et al. (2026) A specialized population of hair afferents dedicated to transmitting mechanical itch. *Neuron*.

Nishinaka T, et al. (2026) IL-4/IL-13 signaling suppresses ABT-263-induced apoptosis in senescent human fibroblasts. *International immunopharmacology*, 184, 116939.

Fischer JS, et al. (2026) A conserved mechanism of membrane fusion in nuclear pore complex assembly. *Molecular cell*.

Chen JJ, et al. (2026) Basal body docking failure triggers centriole clustering and elimination in mammalian spermatogenesis. *EMBO reports*.

Te Marvelde MR, et al. (2026) Assessment of SARS-CoV-2 immune escape using antigenic cartography combined with experimental challenge studies. *NPJ vaccines*, 11(1), 25.

Lassen M, et al. (2026) HR3/ROR γ -mediated cholesterol sensing regulates TOR signaling. *Nature communications*, 17(1).

Sinko DS, et al. (2026) Microporous Immune-Isolating Capsule with Improved Diffusion for Restored Dynamic Bidirectional Hormone Signaling in a Murine Model of Premature Ovarian Insufficiency. *bioRxiv : the preprint server for biology*.

Zheng P, et al. (2026) Evolution of a chloroplast subcompartment housing plastid DNA inside a cup-like starch granule. *Current biology : CB*, 36(8), 2167.

Rivera M, et al. (2026) APOBEC3B enhances the efficacy of PARP inhibitors in elimination of ovarian cancer stem cell. *Scientific reports*, 16(1), 5194.

Qiu F, et al. (2026) Maternal exercise improves vascular function in hypertensive offspring via A-kinase anchoring protein 150 gene (Akap5) epigenetic modifications. *British journal of pharmacology*, 183(10), 2592.

Gerlach I, et al. (2026) Act1LR inhibition improves motor outcome and preserves muscle fibers after experimental autoimmune neuritis. *Acta neuropathologica communications*, 14(1).

Kozono Y, et al. (2026) EC-isHCR: A rapid method for in situ hybridization chain reaction in diverse animal samples. *Methods (San Diego, Calif.)*, 254, 59.

Cheong HSJ, et al. (2026) Organization of circuits linking descending input to motor output in the *Drosophila* Male Adult Nerve Cord connectome. *eLife*, 13.

Mossi Albiach A, et al. (2026) Tissue-reactive zonation of mesenchymal-like phenotypes in human glioblastoma. *Med (New York, N.Y.)*, 101228.

He T, et al. (2026) α -Ketoglutarate couples cellular metabolism to developmental growth via hydroxylation-dependent degradation of Yorkie. *Cell reports*, 45(2), 116922.