

# Resource Summary Report

Generated by [nidm-terms](#) on Jul 29, 2026

## SunCollect MALDI Sprayer

RRID:SCR\_023907

Type: Tool

### Proper Citation

SunCollect MALDI Sprayer (RRID:SCR\_023907)

### Resource Information

**URL:** <https://www.sunchrom.de/en/products/maldi-ms-sample-preparation-2/suncollect-maldi-sprayer.html>

**Proper Citation:** SunCollect MALDI Sprayer (RRID:SCR\_023907)

**Description:** For tissue preparation for Imaging MALDI. The sprayer can handle fresh-frozen or FFPE tissue sections. Enzymes such as trypsin can be applied to the tissue for in situ digestion. Sprayer offers the following key features: Automated delivery of matrix, enzyme or derivatization solutions; 6 port PEEK valve with low dead-volume to minimise liquid usage; Automatic self-cleaning function at end of matrix application; 100 % biocompatible; Compact unit 36 x 49 x 40 cm (w,d,h);Preprogrammed spray protocols;Up to 8 sample slides can be sprayed simultaneously;Easy to use software.

**Synonyms:** SunCollect Sprayer

**Resource Type:** instrument resource

**Keywords:** SunChrom Wissenschaftliche Ger??te GmbH, SunChrom, Tissue Preparation,

**Availability:** Restricted

**Resource Name:** SunCollect MALDI Sprayer

**Resource ID:** SCR\_023907

**Record Creation Time:** 20230802T050218+0000

**Record Last Update:** 20260725T191034+0000

### Ratings and Alerts

No rating or validation information has been found for SunCollect MALDI Sprayer.

No alerts have been found for SunCollect MALDI Sprayer.

---

## Data and Source Information

**Source:** [SciCrunch Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.