

Resource Summary Report

Generated by [nidm-terms](#) on Jul 30, 2026

Children's Hospital of Philadelphia CHOP CAG Genotyping Core Facility

RRID:SCR_028400

Type: Tool

Proper Citation

Children's Hospital of Philadelphia CHOP CAG Genotyping Core Facility
(RRID:SCR_028400)

Resource Information

URL: <https://www.research.chop.edu/cag-genotyping-core>

Proper Citation: Children's Hospital of Philadelphia CHOP CAG Genotyping Core Facility
(RRID:SCR_028400)

Description: Core is equipped with liquid handling technology and experienced staff members. Offers both standard and custom Illumina array project design. With the capacity to process >2300 samples/week and generate >480 million genotypes per day, the genotyping core fully integrates with our Biorepository and Bioinformatics group, offering the collaborator the opportunity for a seamless experience from DNA extraction to data release.

Synonyms: CAG Genotyping Core

Resource Type: access service resource, service resource, core facility

Keywords: ABRF, liquid handling technology, standard and custom Illumina array project design,

Resource Name: Children's Hospital of Philadelphia CHOP CAG Genotyping Core Facility

Resource ID: SCR_028400

Alternate IDs: ABRF_5926

Alternate URLs: <https://coremarketplace.org/?FacilityID=5926&citation=1>

Record Creation Time: 20260508T053733+0000

Record Last Update: 20260729T044713+0000

Ratings and Alerts

No rating or validation information has been found for Children's Hospital of Philadelphia CHOP CAG Genotyping Core Facility.

No alerts have been found for Children's Hospital of Philadelphia CHOP CAG Genotyping Core Facility.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.