

Resource Summary Report

Generated by PRECISE-TBI on Aug 11, 2026

DYKDDDDK-Tag (3B9) Mouse mAb

RRID:AB_2713960

Type: Antibody

Proper Citation

Abmart Cat# M20008, RRID:AB_2713960

Antibody Information

URL: http://antibodyregistry.org/AB_2713960

Proper Citation: Abmart Cat# M20008, RRID:AB_2713960

Target Antigen: DYKDDDDK peptide

Host Organism: mouse

Clonality: monoclonal

Comments: Image validation available for WB, ICC/IF in MDS.

Antibody Name: DYKDDDDK-Tag (3B9) Mouse mAb

Description: This monoclonal targets DYKDDDDK peptide

Target Organism: all

Clone ID: 3B9

Antibody ID: AB_2713960

Vendor: Abmart

Catalog Number: M20008

Record Creation Time: 20250819T233054+0000

Record Last Update: 20250820T053434+0000

Ratings and Alerts

No rating or validation information has been found for DYKDDDDK-Tag (3B9) Mouse mAb.

No alerts have been found for DYKDDDDK-Tag (3B9) Mouse mAb.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 63 mentions in open access literature.

Listed below are recent publications. The full list is available at [PRECISE-TBI](#) .

Xie C, et al. (2026) Lysosome-localized IRTKS condensates promote mTORC1 activity leading to MASLD and HCC. *Cell reports*, 45(2), 116893.

Zhou Z, et al. (2026) The lipid-metabolic enzyme HSD17B12 drives lysosomal degradation of PD-L1 potentiating anti-tumor immunity in a mouse model. *PLoS biology*, 24(1), e3003603.

Chen Q, et al. (2026) Identification of an STING inhibitor targeting the allosteric transmembrane domains. *Cell chemical biology*, 33(5), 623.

Chen YL, et al. (2026) Endometrial stromal cell-derived TMAO sustains decidualization to prevent recurrent spontaneous abortion. *Cell metabolism*.

Song H, et al. (2026) Hijacking ERAD for targeted degradation of transmembrane proteins. *Cell*, 189(6), 1768.

Yan H, et al. (2026) Nanoscale conformational dynamics of human propionyl-CoA carboxylase. *Structure (London, England : 1993)*, 34(1), 62.

Lin J, et al. (2026) AAVLINK: A potent DNA-recombination method for large cargo delivery in gene therapy. *Cell*, 189(3), 969.

Huang R, et al. (2026) TMK-PIN1 drives a short self-organizing circuit for auxin export and signaling in Arabidopsis. *Developmental cell*, 61(1), 73.

Xu F, et al. (2026) Effector-mediated transcriptional rewiring resolves interbacterial conflict through conserved c-di-GMP antagonism. *Current biology : CB*.

Mao Y, et al. (2026) AMPK γ 2 signals amino acid insufficiency to inhibit protein synthesis. *Cell metabolism*, 38(1), 192.

Tang S, et al. (2026) PTEN Loss Promotes PI3K γ Phosphorylation and EPHA2/SRC/p-PI3K γ Y962 Complex Assembly to Drive Tumorigenesis. *Cancer discovery*, 16(3), 521.

Xue D, et al. (2026) Histone methyltransferase SET-18/SMYD2-mediated activation of NADase TIR-1d/SARM1 increases mtROS to promote aging. *iScience*, 29(2), 114649.

Dai P, et al. (2025) Protocol for genetic dissection of class switch recombination using genome-editing tools. *STAR protocols*, 6(2), 103882.

Zheng L, et al. (2025) Phosphorylation-dependent suppression of SERRATE balances miRNA metabolism in planta. *Cell reports*, 44(10), 116367.

Dai P, et al. (2025) Transcription-coupled AID deamination damage depends on ELOF1-associated RNA polymerase II. *Molecular cell*, 85(7), 1280.

Guo S, et al. (2025) TRAF6 regulates ubiquitination-independent TDP-43 condensation and related neurodegeneration. *Molecular psychiatry*.

Tang R, et al. (2025) Stromal Stiffness-Regulated IGF2BP2 in Pancreatic Cancer Drives Immune Evasion via Sphingomyelin Metabolism. *Gastroenterology*, 169(4), 615.

Ding N, et al. (2025) Heat-shock chaperone HSPB1 mitigates poly-glycine-induced neurodegeneration via restoration of autophagic flux. *Autophagy*, 21(6), 1298.

Zhou K, et al. (2025) Antagonistic systemin receptors integrate the activation and attenuation of systemic wound signaling in tomato. *Developmental cell*, 60(4), 535.

Di R, et al. (2025) A mechano-resistance mechanism in skin adapts to terrestrial locomotion. *Cell*.