

Resource Summary Report

Generated by PRECISE-TBI on Jul 31, 2026

CHOP Pathology Core Laboratories

RRID:SCR_009729

Type: Tool

Proper Citation

CHOP Pathology Core Laboratories (RRID:SCR_009729)

Resource Information

URL: <http://eagle-i.itmat.upenn.edu/i/0000013b-9b4b-434b-83a0-df0880000000>

Proper Citation: CHOP Pathology Core Laboratories (RRID:SCR_009729)

Description: Core facility that provides the following services: Tissue microarray construction service, Unstained slides preparation service, H&E staining (on unstained slide) service, Special stain (from slide) service, TUNEL labeling service, IHC stain from slide service, In situ hybridization RNAscope service, New antibody workup service, Decalcification service, Laser capture service, Histopathology-related techniques training, Histopathology technical consultation service, Basic imaging workstation access, Laser capture microdissection microscope access. The Pathology Core Laboratory at the Children's Hospital of Philadelphia Research Institute provides basic histopathology, research immunohistochemistry, tissue microarray, and laser capture microdissection services to researchers at Children's Hospital and within the surrounding academic community. We are located on the 7th Floor of the Leonard and Madlyn Abramson Pediatric Research Center in room 706. The Pathology Core Laboratory unites three core components in a single core facility: histopathology, tissue microarray and laser capture microdissection. The core offers a full range of histopathology services including tissue processing, embedding, and cutting, for both paraffin and frozen tissue. We also perform most standard stains as well as immunohistochemistry, antibody workup, fluorescence, in situ hybridization and TUNEL. Tissue microarrays can be constructed using a Becher Arrayer. Sophisticated imaging instrumentation is available for virtual microscopy (ScanScope from Aperio) and image analysis (Image ProPlus, Volocity). Specialized software is available to image and analyze tissue microarrays, and to manage and store array data.

Resource Type: core facility, service resource, access service resource

Keywords: tissue microarray assay, tissue sectioning, hematoxylin and eosin stain, staining, tunel labeling, in situ hybridization, immunohistochemistry assay, histological sample preparation, laser capture microdissection, light microscopy, histology, fluorescent microscopy

Resource Name: CHOP Pathology Core Laboratories

Resource ID: SCR_009729

Alternate IDs: nlx_156191

Alternate URLs: <http://www.research.chop.edu/cores/pathcore/>

Record Creation Time: 20220129T080254+0000

Record Last Update: 20260731T042818+0000

Ratings and Alerts

No rating or validation information has been found for CHOP Pathology Core Laboratories.

No alerts have been found for CHOP Pathology Core Laboratories.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We found 1 mentions in open access literature.

Listed below are recent publications. The full list is available at [PRECISE-TBI](#) .

Dubensky SB, et al. (2025) Multimodal analysis defines GNG4 as a distinguishing feature of germinal center-positioned CD4 T follicular helper cells in humans. bioRxiv : the preprint server for biology.