

Resource Summary Report

Generated by [PRECISE-TBI](#) on Aug 1, 2026

Rigaku: X-Stream 2000 Cryo-Freezing System

RRID:SCR_018633

Type: Tool

Proper Citation

Rigaku: X-Stream 2000 Cryo-Freezing System (RRID:SCR_018633)

Resource Information

URL: <https://www.labx.com/item/rigaku-x-stream-2000-cryogenic-crystal-cooler-system/LV40211010>

Proper Citation: Rigaku: X-Stream 2000 Cryo-Freezing System (RRID:SCR_018633)

Description: THIS RESOURCE IS NO LONGER IN SERVICE. Documented on March 11,2024. Self contained cryo freezing system designed cooling crystals with range of negative 80C to 4C. It produces nitrogen, which is then cooled to negative 196C, and then rewarmed to desired temperature.

Resource Type: instrument resource

Keywords: Cryo freezing System, Instrument, Equipment, USEDit, Rigaku, ABRF

Availability: THIS RESOURCE IS NO LONGER IN SERVICE

Resource Name: Rigaku: X-Stream 2000 Cryo-Freezing System

Resource ID: SCR_018633

Alternate IDs: Model_Number_X-Stream_2000

Record Creation Time: 20220129T080341+0000

Record Last Update: 20260725T190851+0000

Ratings and Alerts

No rating or validation information has been found for Rigaku: X-Stream 2000 Cryo-Freezing System.

No alerts have been found for Rigaku: X-Stream 2000 Cryo-Freezing System.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.