

Resource Summary Report

Generated by PRECISE-TBI on Aug 3, 2026

Dusk2Dawn

RRID:SCR_027417

Type: Tool

Proper Citation

Dusk2Dawn (RRID:SCR_027417)

Resource Information

URL: <http://www.nitrc.org/projects/dusk2dawn/>

Proper Citation: Dusk2Dawn (RRID:SCR_027417)

Description: Software plugin for 'EEGLAB' MATLAB toolbox, which allows to whole-night sleep EEG using Artifact Subspace Reconstruction (ASR). Used to implement Artifact Subspace Reconstruction cleaning with whole-night sleep EEG. Automatic cleaning of whole-night sleep electroencephalogram using Artifact Subspace Reconstruction.

Resource Type: software resource, software application

Keywords: whole-night sleep electroencephalogram, automatic cleaning of whole-night sleep electroencephalogram, Artifact Subspace Reconstruction, implement Artifact Subspace Reconstruction cleaning, whole-night sleep EEG,

Availability: Free, Available for download, Freely available,

Resource Name: Dusk2Dawn

Resource ID: SCR_027417

Record Creation Time: 20250912T055219+0000

Record Last Update: 20260801T191421+0000

Ratings and Alerts

No rating or validation information has been found for Dusk2Dawn.

No alerts have been found for Dusk2Dawn.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.