

Resource Summary Report

Generated by [PRECISE-TBI](#) on Aug 6, 2026

Trikinetics DAM FileScan

RRID:SCR_021817

Type: Tool

Proper Citation

Trikinetics DAM FileScan (RRID:SCR_021817)

Resource Information

URL: <https://manualzz.com/doc/6543996/drosophila-activity-monitoring-system-user-s-guide>

Proper Citation: Trikinetics DAM FileScan (RRID:SCR_021817)

Description: Software tool used in Drosophila Activity Monitor system to scan data for errors. Also allows user to extract data from monitor files and bin data in different lengths.

Synonyms: DAM FileScan

Resource Type: software resource

Keywords: scan data for errors, Drosophila Activity Monitor system, extract data

Resource Name: Trikinetics DAM FileScan

Resource ID: SCR_021817

Record Creation Time: 20220129T080357+0000

Record Last Update: 20260801T190718+0000

Ratings and Alerts

No rating or validation information has been found for Trikinetics DAM FileScan.

No alerts have been found for Trikinetics DAM FileScan.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We found 1 mentions in open access literature.

Listed below are recent publications. The full list is available at [PRECISE-TBI](#) .

Iyengar AS, et al. (2022) Under warm ambient conditions, *Drosophila melanogaster* suppresses nighttime activity via the neuropeptide pigment dispersing factor. *Genes, brain, and behavior*, 21(4), e12802.