

Resource Summary Report

Generated by SPARC SAWG on Sep 14, 2026

Select Science: Genomic Solutions OmniGrid 100 Microarrayer

RRID:SCR_018073

Type: Tool

Proper Citation

Select Science: Genomic Solutions OmniGrid 100 Microarrayer (RRID:SCR_018073)

Resource Information

URL: <https://www.selectscience.net/products/omnigrid-100/?prodID=11221>

Proper Citation: Select Science: Genomic Solutions OmniGrid 100 Microarrayer (RRID:SCR_018073)

Description: OmniGrid 100 microarrayer is high capacity multi axis microarrayer that prints biological samples onto solid substrates with speeds that give high throughput array production. Can be integrated with optional Server Arm 3 for fully automated handling of up to 72 input plates.

Synonyms: GeneMachines OmniGrid 100 Microarray Printer, OmniGrid 100 microarrayer

Resource Type: instrument resource

Keywords: ABRF, microarrayer, printer, instrument, equipment, Genomics solutions Ltd.

Resource Name: Select Science: Genomic Solutions OmniGrid 100 Microarrayer

Resource ID: SCR_018073

Alternate IDs: Model_Number_OmniGrid 100

Alternate URLs: <https://www.accela.eu/files/products/97/40domnigridmicro.pdf>

Record Creation Time: 20220129T080338+0000

Record Last Update: 20260912T195857+0000

Ratings and Alerts

No rating or validation information has been found for Select Science: Genomic Solutions OmniGrid 100 Microarrayer.

No alerts have been found for Select Science: Genomic Solutions OmniGrid 100 Microarrayer.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We found 1 mentions in open access literature.

Listed below are recent publications. The full list is available at [SPARC SAWG](#) .

Xiao W, et al. (2026) Metabolite-gated vascular contractility switch: OXGR1 activation mechanism enables agonist therapy for rosacea erythema. Cell, 189(7), 1990.